



There's something missing
at the core of your core lab.



Once in a while, a new product comes along
that changes how people work.

For core labs, the **Fragment Analyzer™** provides unmatched versatility.

- ◆ Analyze samples ranging from 35bp to 40,000bp
- ◆ Perform quality and quantity checks on DNA fragments prior to Next Gen sequencing runs
- ◆ Measure total RNA with RQN number, which assesses total RNA integrity on a scale of 1 to 10 and correlates to legacy data
- ◆ Eliminate bottlenecks from slab gel and lab-on-chip methods
- ◆ Available in 12 or 96 sample formats

**“The Fragment Analyzer is easy to use
and performs many tasks.”**

—Dr. Caroline Allix-Béguec
R&D Manager, Genoscreen, Lille, France



info@aati-us.com www.aati-us.com

Call our lab in Heidelberg to learn more: 49-6221-868058-10

