



High-throughput RNA analysis



RQN provides quality measurement

RNA research is the **new frontier** of discovery in life science. Labs that once analyzed dozens of RNA samples a week may now analyze hundreds or even thousands.

The Fragment Analyzer's **single dilution step** for RNA sample prep, coupled with parallel automated separation, achieves high-throughput. One dilution – in either **12- or 96-sample formats** – and it's ready to go.

RQN means **RNA Quality Number**. Derived from PROSize® analytical software, RQN assesses the quality or **integrity** of total RNA on a scale of **1-to-10**, and allows for **correlation to legacy data**.

Every once in a while, a new product comes along that changes how people work. *Find out for yourself.*



Phone: 515-296-6600 | www.aati-us.com